



REACH DEEP-WELL 384 OR 1536 PLATES

BIOROBOTIX™ 384-WELL EXTENDED LENGTH 50UL

TIPS DESIGNED TO FIT THE BIOMEK FX/NX & MULTIMEK AP384 WORKSTATIONS

- Increased volume and high throughput screening with extended ability to reach into deep-well 384 and 1536 plates
- Pipette tips, trays and tray lids reduce static charges so tips pick up straight and eject perfectly each time
- Tips are lot-tested on the actual workstations for which they are created using 15-step quality control tests
- Utilizes the most of the Beckman workstation and 384-well head
- 384-well Extended Length 50µl tips will be available as sterile, non-sterile and with our patented ART® barrier



Cat. #	Description	Packaging
938-261	384-well 50ul Extended Length Clear Sterile Tips	96 tips/tray, 10 trays/pack, 5 packs/case
938-262	384-well 50ul Extended Length Clear Non-Sterile Tips	96 tips/tray, 10 trays/pack, 5 packs/case
938-021	384-well 40ul Extended Length Clear ART® barrier Tips	96 tips/tray, 10 trays/pack, 5 packs/case

**CONTACT US FOR YOUR FREE
VALIDATION SAMPLES TODAY!**

MBP[®]
Molecular BioProducts
THE INTELLIGENT CHOICE FOR DISCOVERY.

9389 Waples Street, San Diego, CA. 92121 U.S.A
 Within the U.S. call 800.995.2787 or Fax 858.453.4367
 Outside the U.S. call 858.453.7551
 email: info@mbpinc.com www.mbpinc.com

MBP is a registered trademark of Molecular BioProducts, Inc. ART is a registered trademark of MBP and is protected by U.S. Patent No. 5,156,811.

08082407

