

FOR IMMEDIATE RELEASE

Contact: Mitch Becker
Phone: (858)453-7551 ext.152

Molecular BioProducts® Introduces Genomic Tips for Automated Workstations

San Diego, CA – Molecular BioProducts, Inc., a market-leading manufacturer of pipet tips for automated workstations introduces new genomic tips for Multimek Type II, QIAGEN/Rosys, Tecan and MultiPROBE II HT automated workstations. These products compliment our existing genomic tips offered for Beckman Biomek FX, Multimek AP96, and Biomek 1000/2000 workstations. This new style of tip is the same high quality MBP® BioRobotix™ brand tip but with a wide-bore orifice (genomic). These genomic tips are perfect for pipetting viscous tissue samples such as blood as well as for preventing damage to genomic DNA during sample preparation and processing. All tips come in a variety of sizes as well as being available as sterile, non-sterile and with the patented ART® barrier for guaranteed protection against crossover contamination.

The MBP® BioRobotix™ brand tips is the most comprehensive product line available for automated workstations with tips to fit most major workstations. For more information, contact Molecular BioProducts at (800)955-2787 or (858)453-7551. Additionally, we may be contacted by sending an e-mail to info@mbpinc.com or visiting our website at www.mbpinc.com.