

MBP[®] Molecular BioProducts

FOR IMMEDIATE RELEASE

Contact: Mitch Becker
Phone: (858) 453-7551 ext.152

Molecular BioProducts[®] Introduces the 1000 REACH Pipet Tip

San Diego, CA – February 2002, Molecular BioProducts, Inc., a market-leading manufacturer of pipet tips for molecular biology introduces its newest pipet tip, the 1000 Reach. Expanding on the success of one of its most popular pipet tips, the 10ul Reach, MBP has created a truly unique 1ml pipet tip that is double the length of many standard 1ml pipet tips.

This extended length pipet tip is nearly five inches long, and has been designed to reduce the chances of sample and pipettor contamination by reaching to the bottom of deep storage containers such as 15ml centrifuge tubes without the pipettor shaft entering the sample vessel.

This tip also features a MicroPoint design to reduce sample hold up and SoftFit technology to ensure a secure, easy fit on all 1ml pipettors. Molded from virgin polypropylene, this tip is available with and without the patented ART self-sealing barrier and is certified to be RNase, DNase, Nucleic Acid and Pyrogen Free.

For more information, contact Molecular BioProducts at (800)955-2787 or (858)453-7551. Additionally, we may be contacted by sending an e-mail to info@mbpinc.com or visiting our website at www.mbpinc.com.



ART 1000 Reach Pipet Tip