



**FOR IMMEDIATE RELEASE**

Contact: Kelly D. Halley  
Phone: (858) 453-7551 ext.159

## **Molecular BioProducts® Introduces Our New BioRobotix Biomek FX Span-8 Tips**

San Diego, CA, October 2002 – Molecular BioProducts, Inc., a market-leading manufacturer of pipet tips for automated workstations introduces our newest product: liquid sensing (conductive) 96-well tips for Beckman Biomek FX workstations with a **Span-8 pod**. We have taken our years of success in molding top quality Tecan and Qiagen liquid sensing tips and created the perfect design for the Span-8 pod. MBP is the first tip manufacturer other than Beckman to offer these tips and the only manufacturer who will offer them in sterile, non-sterile, and filtered versions (with our patented ART® filter).

These new tips, like all lots of our BioRobotix™ brand tips, are quality control tested on the actual workstation for which they were created. These 15 quality control tests include parameters such as Total Indicator Runout (TIR) or tip straightness and Percentage Coefficient of Variance (%CV). Molecular BioProducts is the only tip manufacturer that carries out these tests on the actual workstations before release.

MBP® BioRobotix™ brand is the most comprehensive product line of tips available for automated workstations on the market today. For more information or free samples of this product, contact Molecular BioProducts at (800)955-2787 or (858)453-7551. Additionally, we may be contacted by sending an e-mail to [info@mbpinc.com](mailto:info@mbpinc.com) or visiting our website at [www.mbpinc.com](http://www.mbpinc.com).

