

Contact Bill LeMaster  
Phone: (858) 453-7551 ext.157

## **Molecular BioProducts Announces ART® REACH Pipet Tips**

San Diego, CA, October 2003, Molecular BioProducts, (MBP), leading innovator in liquid handling technology offer a unique family of pipet tips. The one-of-a-kind REACH pipet tips adds a layer of security in protecting samples and pipettors by preventing the pipettors shaft from touching the inside of the sample vessel, thus virtually eliminating the chance of carry over contamination and the need to clean or replace pipettors. Assays remain pure, saving time and effort.

MBP's 1000 REACH is nearly 5 inches long, twice the length of regular 1ml pipet tip. It is designed to reach the bottom of 15ml & 50ml centrifuge tubes and other deep sample vessels. Also available with the ART® self-sealing barrier, REACH pipet tips are offered in 10ul, 100ul, 200ul, and 1000ul volumes.

For more information or to request a free sample, log on to Molecular BioProducts' web site at [www.mbpinc.com](http://www.mbpinc.com), send an email to [info@mbpinc.com](mailto:info@mbpinc.com), or contact us at (800) 955-2787 or (858) 453-7551.



**ART 1000 Reach, 1000ul Extended Length, Aerosol Resistant Tip (ART)**

ART is a registered trademark of MBP and is protected by U.S. Patent Number 5,156,811;  
MBP is a registered trademark of Molecular BioProducts, Inc., San Diego, CA.